

Wizfill

Wizfill with Lidocaine 1.1mL



- ✓ **23G Cannula**
(Only Wizfill Lidocaine No.2, No.3)
- ✓ **27G Needle**
- ✓ **0.3% Lidocaine**
- ✓ **24mg/mL HA Contents**
- ✓ **Monophasic Filler**

Wizfill fills up all the aged signals by proper volumes and natural contours.

- + Monophasic Gel¹⁾
- + 0.3% Lidocaine

Wizfill

Cross Linked Hyaluronic acid Dermal Filler, Wizfill



Wizfill No.1 with Lidocaine 1.1mL

Ingredient

Lidocaine hydrochloride.....3mg/mL
Cross-linked hyaluronic acid...24mg/mL
Phosphate buffered saline, pH 7.....q.s.

Volume

1.1mL× 1syr

Needle

27G × 2EA

Injection Depth

Superficial Dermis

Wizfill No.2 with Lidocaine 1.1mL

Ingredient

Lidocaine hydrochloride.....3mg/mL
Cross-linked hyaluronic acid...24mg/mL
Phosphate buffered saline, pH 7.....q.s.

Volume

1.1mL× 1syr

Needle / Cannula

27G × 1EA / 23G × 1EA

Injection Depth

Superficial to Mid Dermis

Wizfill No.3 with Lidocaine 1.1mL

Ingredient

Lidocaine hydrochloride.....3mg/mL
Cross-linked hyaluronic acid...24mg/mL
Phosphate buffered saline, pH 7.....q.s.

Volume

1.1mL× 1syr

Needle / Cannula

27G × 1EA / 23G × 1EA

Injection Depth

Deep Dermis to Subcutaneous

Wizfill

Specific Treatment Area



Forehead

Wizfill No.2 with Lidocaine 1.1mL
Wizfill No.3 with Lidocaine 1.1mL

Nose

Wizfill No.3 with Lidocaine 1.1mL

Cheek

Wizfill No.1 with Lidocaine 1.1mL
Wizfill No.2 with Lidocaine 1.1mL

Nasolabial folds

Wizfill No.1 with Lidocaine 1.1mL
Wizfill No.2 with Lidocaine 1.1mL

Lip

Wizfill No.1 with Lidocaine 1.1mL

Chin

Wizfill No.3 with Lidocaine 1.1mL

Cross Linked Hyaluronic acid Dermal Filler, Wizfill

Wizfill

** In addition to the models recommended for each injection site, other models are also available.*



WIZMEDI

WIZMEDI CO., LTD.

#216, B/D.No.G12, 280, Daehak-ro, Gyeongsan-si, Gyeongsangbuk-do, 38541, Republic of Korea

www.wizmedi.kr | **MADE IN KOREA**